

Date of Document: 18/3/2025
Clinic Name: Stride Podiatry
Clinic Address: 66 McLennan St Albion 4010 QLD
Clinic Phone: 0.46919798
Clinic Contact: Cristina

Sterile instrument order:

| Item | Quantity | Sterilisation Date | Tracking number | Expiry Date |
|----------|----------|--------------------|-----------------|-------------|
| Galipots | 5 | 17/3/2025 | 250317:1 | 16/3/2026 |
| | | | | |
| | | | | |
| | | | | |
| | | | | |
| | | | | |
| | | | | |

Total items: 5 **Batch** 250317:1 **Sets** 20
Clean boxes Supplied: 1
Dirty Boxes supplied: 1

Sterile compliant on autoclave
 Sterile compliant on pack
 Sterile compliant on class 6 emulating indicator
 Sterile compliant on indicator label

YES / NO
 YES / NO
 YES / NO
 YES / NO

Certified By: Chris H

Signed: 

CYCLE DATE 17/03/25 AUTOCLAVE 1
~~BROWN-STERILIZED~~
 250317-1
 BATCH NUMBER

Emulating Indicator
 134 °C - 3,5 min. / 121 °C - 15 min.

REF: 3FSKS645618
 EN-ISO 11140-1
 Type 6 Indicator
 LOT: 1318267
 2023/05 2026/05
 Manufactured in The Netherlands.
 www.interster.nl
INTERSTER

STEAM
 Minimum reference colour
 Sterilizing conditions turn indicator from light blue to dark.

Assembled by: _____
 Date: _____

Notes: